

Technical Specification

Product Name	ATBW-110
Category	Anesthesia Trolley
Measurements and Weight	
Dimensions (W x D x H) [mm]	110 x 87 x 59,5
Net Weight [kg]	49
Max. Total Loading [kg]	200
Equipment	
Storage	
No. of Drawer	5
Measurement of the Drawers	1. and 2. Drawers : 570x440x65mm 3. and 4. Drawer: 570x440x135mm 5. Drawer:570x440x220mm
Mobility aids	
Castors	4 single castors (Ø 125 mm), two with brakes and one with directional lock
Wall Bumper	Standard
Materials	
Pillar	ABS Plastic
Side and Back Boards	ABS Plastic
Worktop	ABS Plastic
Drawer	ABS Plastic
Drawer handle	ABS Plastic
Wall Bumper	ABS Plastic
Color	Coince of 4 colors, combined freely
Accessories	
Medicine Boxes	Standard
Central Lock System	Standard
Extendable Table	Standard
Soft-Closing Drawer Runners	Standard
Waste Bin	Standard
Drawer Divider, transparent	Standard
Defibrillator Shelf	Optional
CPR Board	Optional
Oxygen Tank Holder	Optional
IV Hang (2 pole)	Optional
Power Socket	Optional
Disposable Needle Holder	Optional