

TECHNICAL DOCUMENTATION

Product Name	Emergency Cart
Model	ECRW-110 / ECRG-110
Rev. No.	00
Page No	1 / 1



TECHNICAL SPECIFICATIONS

Characteristic	Explanation
Body Material	Durable ABS plastic
Color Options	Red-white or red-grey
External Dimensions	1100 × 870 × 600 mm
Carrying Capacity	>100 kg
Wheel	4 units, 125 mm diameter, 2 with brake
Lock System	Single-use central lock (standard) + optional reusable lock
Number of Drawers	Total 5 pieces: 2 small, 2 medium, 1 large
Drawer Structure (Optional)	On request, the number of drawers can be changed during the production phase
Drawer Dimensions	2 x Small: 61×48×8 cm, 2 x Medium: 61×48×16 cm, 1 x Large: 61×48×25 cm
Interior Partition Options	Small: 24/48 bays, Medium: 12 bays, Large: 24 bays
Label Holder	5 transparent drawer labels
Openable Surface	Drug preparation table
Handles	Ergonomic guiding handles
Cleaning Resistance	Resistant to disinfectants and UV, antibacterial surface

STANDARD ACCESSORIES

Accessory	Explanation
Waste Bin	10 lt ABS
Material Basket	ABS, 3-compartment
Defibrillator Stand	Metal
Locks	10 disposable locks
Label Holder	5 pcs.
Serum Strap	2-pack serum strap
CPR Board	ABS material
O2 Tube Mounting Kit	ABS or metal (for 1–3–5 lt)

OPTIONAL ACCESSORIES

Accessory	Explanation
Double Waste Bin	ABS, 2-compartment
ABS Defibrillator Stand	Alternative structure
Reusable Lock	Mechanism system
4-Piece Serum Strap	Large carrying capacity
Power Outlet	3-piece socket
Mobile Aspirator Rack	Metal material
Needle Waste Bin	For injection safety

WARRANTY AND SERVICE

Title	Explanation
-------	-------------

TECHNICAL DOCUMENTATION	Product Name	Emergency Cart
	Model	ECRW-110 / ECRG-110
	Rev. No.	00
	Page No	2 / 1

Warranty Period	Minimum 1 year
Spare Parts Supply	Minimum 5 years
DELIVERY AND INSTALLATION	
Title	Explanation
Delivery Time	30 days from the date of order
Packaging	Foam or carton box, impact proof
Installation Status	Installed delivery or easy installation is in structure